



i<sup>™</sup> seri  
Non-i i ilators

**Patient name**  
(please print)

Lumis 150 ST-A iVAPS therapy settings (pressure range of 4-30cm/H<sub>2</sub>O)



**Mode**  iVAPS

**Patients Height** \_\_\_\_\_ cm **Target respiratory rate** \_\_\_\_\_ bpm (8-30 bpm) **Target Va** \_\_\_\_\_ L/min

**EPAP** \_\_\_\_\_ cm/H<sub>2</sub>O - Please complete or NOT using Auto EPAP

**Auto EPAP** **Min EPAP** \_\_\_\_\_ cm/H<sub>2</sub>O **Max EPAP** \_\_\_\_\_ cm/H<sub>2</sub>O - Please complete if using Auto EPAP

**Min PS** \_\_\_\_\_ cm/H<sub>2</sub>O **Max PS** \_\_\_\_\_ cm/H<sub>2</sub>O **Ti Max** \_\_\_\_\_ sec **Ti Min** \_\_\_\_\_ sec

**Rise Time** \_\_\_\_\_ ms (100-900ms)

**Trigger sensitivity**  Very low  Low  Medium  High  Very high  
Less sensitivie/Harder to trigger Default More sensitivie/Easier to trigger

**Cycle sensitivity**  Very low  Low  Medium  High  Very high  
Less sensitivie/Harder to cycle Default More sensitivie/Easier to cycle

**Mask**  Full face  Nasal  Pillows  Other - Please specify: \_\_\_\_\_

**Ramp**  Off  On **Specified time** \_\_\_\_\_ min **Start EPAP** \_\_\_\_\_ cm/H<sub>2</sub>O Default is 4cm/H<sub>2</sub>O  
Note: Ramp can be set between 5-45mins in increments of 5mins

**Does the patient use supplemental oxygen with their device**  No  Yes \_\_\_\_\_ L/min  
Note: the maximum supplemental oxygen flow into this device in iVAPS mode is 4 L/min

**Alarms**

**Low MV** \_\_\_\_\_ L/min  Off **High leak**  Off  On **Non vented mask**  Off  On

**Low SpO<sub>2</sub>**  Off  On \_\_\_\_\_ % (if SpO<sub>2</sub> monitor is attached) **Apnea alarm** \_\_\_\_\_ Sec  Off

**Alarm volume**  Low  Medium  High