

# Removable Rigid Dressings / Rigid Casts

## Nursing Information

Patient's sticker

This patient has been fitted with a removable rigid dressing / rigid cast.


The cast was fitted on (date):

The dressings present at the time of casting were:

The number of socks applied at the time of fitting:

The cast is held on using a suspension sock.

- Please keep the profile / thickness of the wound dressing consistent as it will affect the fit of the cast.
- Do not remove the rigid dressing for more than 10 minutes at a time.
- Wool socks should be added as the volume of the residual limb decreases.
- When the patient requires 3 or more wool socks or if there are problems fitting the rigid cast, it is time to have a new cast fitted.



- TUBIGRIP SOCK
- STUMP SOCKS
- RIGID CAST