

DreamStation BiPAP AVAPS

Page 1 of 2

PHILIPS

Sleep and Respiratory Care

65 Epping Road,
North Ryde NSW 2113

1800 830 517

clinical-src@philips.com

www.philips.com.au/respironics

Prescription form for pressure modes –
no AVAPS, no AAM

Name

Date of birth

Address

Phone



DreamStation BiPAP AVAPS is intended to provide mask-applied non-invasive ventilator support to Obstructive Sleep Apnoea (OSA) and Respiratory Insufficiency patients weighing over 18kg. This device may be used in the hospital or the home.

Modes

CPAP

Spontaneous Timed (S/T)

Spontaneous (S)

Timed (T)

Pressure-Control (PC)

Settings

CPAP Pressure	cmH ₂ O	Ramp	On	Off
<i>CPAP Mode Only</i>				
4 – 20cmH ₂ O		Ramp time (5 – 45 min)		
(1cmH ₂ O increments)		Min (in increments of 5 mins)		
		Ramp start (4 – 20cmH ₂ O)		
		cmH ₂ O (in increments of 0.5cmH ₂ O)		
IPAP Pressure cmH ₂ O	cmH ₂ O			
<i>Min 4 – Max 30cmH₂O</i>				
<i>(0.5cmH₂O increments)</i>				
EPAP Pressure	cmH ₂ O			
<i>Min 4 – Max 25cmH₂O</i>				
<i>(0.5cmH₂O increments)</i>				
Bi-Flex Setting	1,2,3			
<i>(S mode only – rise time disabled)</i>				
Bi-Flex Lock	On	Off		
Back Up Rate	bpm			
<i>S/T, T & PC Modes Only</i>				
<i>Off, 1 – 30 BPM</i>				
<i>When in Timed mode, the minimum setting is 4 BPM</i>				
Inspiratory Time	s			
<i>(S/T, T, PC modes only) 0.5 to 3 seconds in 0.1 increments</i>				
Rise Lock	On	Off		
Rise Time				
<i>In S mode, Bi-Flex must be OFF</i>				
<i>1 – 6 in increments of 1</i>				

Humidification

Off

On – default settings

Humidifier 3, Tube Temperature 3

On – non-default

Please complete all settings below

(Note: if settings incomplete default settings will be applied)

Humidification Adaptive Fixed

Humidifier 0 (off), 1-5 (1 increment)

Tube temperature 0 (off), 1-5 (1 increment)

clinical-src@philips.com

Page 1 of 2

Interface	Supplementary Oxygen	No	Yes	Flow Rate	L/min Max 15 L/min
Circuit Type	Optional – Standard 15mm heated tube provided				

Alarms (Mandatory - Select On or Off for each alarm type)

Alarm Type	Off	On	Alarm Setting	Alarm Parameter
Patient Disconnect Alarm			sec	Off, 15 or 60 sec
Apnoea Alarm			sec	Off, 10, 20, 30, 40 sec
Low Minute Vent Alarm			LPM	Off, 1 – 99 LPM in 1.0 LPM increments
Low Tidal Volume Alarm				AVAPS only. Will sound if IPAP Max is reached and Vt is not achieved.

Connected Care

Add Prescriber to Care Orchestrator Care Team Cellular Modem Use DreamMapper Mobile App

Additional Information

Physician Name	Email
Physician Signature	Phone
Date	