

prisma30ST Prescription Form



Fill out one prescription per form.

Lock form at the end once completed.

Please note: Some web browsers may render form fields differently. For best results, download script and open with Adobe Reader.

Patient Information

Patient Name	Unit/House No.
Date of Birth	Street
Patient ID/MRN	Suburb
	State
Date	Post Code

Device Settings

Language Expert	English	Other	Language
Language Patient	English	Other	Language
Brightness	(1-3)		
Energy saving	Off	On	
Patient info menu	All	Usage	
Key tone volume	(0-3)		
Alarm Volume	(1-3)		

Comfort Settings

prisma AQUA Humidifier	No	Yes	
Smart Aqua Control	On	Off	
Humidifier Level	(1-7)		
Humidifier Level max	(1-7)		
Mask test pressure	(4-20 cmH2O)		

Accessories

Supplemental Oxygen (via circuit connector)	No	Yes	(up to 15 l/min)
Bacterial Filter	No	Yes	
Tube Type	19-22mm	15mm	
Tube Type Lock	Off	On	
Air Filter Reminder	Off	On	

prisma30ST Prescription Form

	Setting Range	Increment	CPAP	APAP	S	ST	autoST	aPCV	T
Select mode									
CPAP/P min	4 - 20 cmH2O	0.5 cmH2O							
P max	4 - 20 cmH2O	0.5 cmH2O							
APAP	std, dyn	-		std dyn					
IPAP/IPAP Min*	4 - 30 cmH2O	0.5 cmH2O							
EPAP/EEPAP*	4 - 30 cmH2O	0.5 cmH2O							
EPAP max/EEPAP max**	4 - 25 cmH2O	0.5 cmH2O							
EPAP min/EEPAP min**	4 - 25 cmH2O	0.5 cmH2O							
PDIFFF**/PDIFFF min*	0 - 26 cmH2O	0.5 cmH2O							
Target Volume									
Off									
On									
Target Volume Additional Settings (Only complete if Target Volume is selected On)									
Target Volume	300 - 2000ml	10ml							
Pressure Modification	1 - Slow 2 - Medium 3 - Fast	1 - 0.5/8 breaths 2 - 1.0/5 breaths 3 - 1.5/breath							
IPAP max*	4 - 30 cmH2O	0.5 cmH2O							
PDIFFF max*/**	0 - 26 cmH2O	0.5 cmH2O			**	**	**/**	**	
softPAP									
Off									
softPAP 1									
softPAP 2									
softPAP 3									
softPAP Lock			Off	On	Off	On			
Bi softPAP									
Off									
Bi soft 1									
Bi soft 2									
TriLevel									

*Target Volume Setting
**TriLevel Setting

prisma30ST Prescription Form

	Setting Range	Increment	CPAP	APAP	S	ST	autoST	aPCV	T
Settings 2 - Back up rate and Inspiratory Time									
autoF	On, Off	-					On Off		
F min/F (Back-up Rate)	Off, 5 - 35 bpm aPCV, T: 5 - 35 bpm	1 bpm							
Ti Timed	Auto or 0.3 - 4s	0.1s							
Ti Max	0.3 - 4s	0.1s							
Ti Min	0.3 - 4s	0.1s							
Ti	0.3 - 4s	0.1s							
Settings 3 - Ramp and Trigger									
Ramp IN	0, 1, 2, 3	-							
Ramp EX	1, 2, 3	-							
Trigger IN	1, 2, 3 or Auto	-							
Trigger EX	1, 2, 3 or Auto	-							
Comfort - autoSTART									
Off									
On									
autoSTOP									
Off					Locked Off	Locked Off	Locked Off	Locked Off	Locked Off
On									
softSTART									
Off									
On									
softSTART Settings (Only complete if softSTART is selected On)									
softSTART max	5 - 45 mins	5 mins							
softSTART pressure min	4 - 20 cmH2O	0.5 cmH2O							
softSTART PDIFF	Off, On	Off, On			Off On	Off On	Off On	Off On	
softSTOP									
Off									
On									
softSTOP Settings (Only complete if softSTOP is selected On)									
softSTOP max	5 - 45 mins	5 mins							

prisma30ST Prescription Form

			Alert Off/On	CPAP	APAP	S	ST	autoST	aPCV	T
	Setting Range	Increment								
Alerts	(If alert is selected On, provide parameter where applicable in the text box)									
Disconnection	Off, On	-	Off On							
Severe Leak	Off, On	-	Off On							
Leakage Lock Alert	Off, On	-	Off On							
Low VT	Off, 300ml - 2000ml	10ml	Off On							
Low MV	Off, 1 - 25l/min	1 l/min	Off On							
Long Apnea	Off, 10 - 60s	10s	Off On							

prisma30ST Prescription Form

Additional Information

Prescribers Details

Name Signature

Health Provider

provider No.

Phone Number

E-mail

CHECK BOX TO LOCK THE FORM (cannot be unticked)

Technical Setup Sign Off (Internal Use only)

Device Type Humidifier Serial Number

Device Serial Number Warranty End Date

Technician 1

Name

Phone Number

Signature

Date

Time

Technician 2

Name

Phone Number

Signature

Date

Time

Additional Information